



The Medical Power of Light

Facial Aesthetics with the LightWalker[®] Dental Laser

Now available for your dental practice

- Safe and patient-friendly procedures
- Minimally invasive and tissue-selective effect
- Natural results
- New treatment possibilities
- Gentle, fast and convenient
- Generates additional income

Fotona
choose perfection

SmoothEye®

A gentle Er:YAG laser treatment for tightening and wrinkle reduction in the periorcular region.

Before

After



Video



Courtesy of Jolanda Grayling
Laser source: Non-Ablative Er:YAG

Vascular lesions

Instant results can be achieved by treating vascular lesions with Nd:YAG.

Before

After



Video



Courtesy of Robin Sult

Fibroma removal

Precise and fast removal of benign lesions with Er:YAG.

Before

After



Video



Courtesy of Dr. Dovsak

Hair removal


Safe and effective, with long-lasting results for all skin types using Nd:YAG.

Before

After



Video



Courtesy of Robin Sult



Fotona3D™

Non-invasive 3-step laser facelift with a synergistic combination of the Er:YAG and Nd:YAG wavelengths.

Before



Video



After



Courtesy of Dr. S. Koceva

- 1st step: SmoothLiftin® (intraoral mucosa tightening)
- 2nd step: Long-pulse Nd:YAG nonablative skin rejuvenation
- 3rd step: Er:YAG gentle skin resurfacing

LipLase®

Fast and patient-friendly lip enhancement treatment with Fotona SMOOTH® mode.

Before



Video



After



Courtesy of Dr. Harvey Shiffman

Skin resurfacing

Multiple Er:YAG treatment modalities for enabling the clinical outcomes your patients desire.

Before



Video



After



Courtesy of Dr. L. Volovec




"To be a dentist in the 21st century you have to get out of the mouth."


— Dr. Linhlan Nguyen, BDS, MSc

LightWalker® lasers are designed for all forward-thinking professionals who wish to upgrade their daily dental experience with new treatment possibilities that only the latest laser technology can offer.

Nd:YAG Lasers – for Treatment of Deep-Lying Structures



Er:YAG Lasers – for Treatment of Superficial Lesions




Fulfill your patients' expectations with a wide range of treatment options, from cold and very superficial skin resurfacing to deep and hot interaction with the tissue. Additionally, you can offer non-invasive skin and mucosa tightening with Fotona SMOOTH® mode.

Advantages of Fotona SMOOTH® Mode

Fotona's revolutionary non-ablative SMOOTH® mode treats the skin and mucosa in a smooth, almost "feather-like" non-ablative manner, without bleeding and with precisely controlled temperature deposition, making it ideal for skin tightening. The intense, controlled surface tissue heating stimulates collagen remodeling and initiates neocollagenesis. The end result is an overall improvement of wrinkles, laxity and elasticity in the treatment areas.

SMOOTH mode pulse



Getting started in facial aesthetics with LightWalker®

Training in facial aesthetics is provided through the Laser and Health Academy (www.laserandhealth.com) under the guidance of highly skilled experts in medical laser technology, as part of a training workshop covering Fotona's minimally invasive aesthetic treatments. The extensive workshop, where participants engage in live demonstrations and gain an in-depth understanding of laser-tissue interaction, will provide the needed insight into the fundamentals of non-invasive aesthetic treatments and other procedures that can be performed with Fotona laser systems.

To learn more about Facial Aesthetics and what LightWalker® can do for your practice, contact Fotona at info@fotona.com today.